



MARCA LAB

ACUPUNCTURE ▶

SEIRIN®
Shōnishin
Moxa

Acupuncture Models
Cupping Sets
Massage & Wellness

3B LASER
Taping

Wall Charts
Massage Tables

ACUPUNCTURE

Dear customer,

To provide your patients with the highest standard of treatment, you rely on high quality equipment for your practice. On the following pages, you will find a selection of products specifically arranged for the needs of successful therapy, that will save you time and optimise your treatments.

The new SEIRIN® M-Type uses cutting-edge needle technology with its guide tube and a longer handle for a better grip. Please see page 11 for more information. SEIRIN® stands for high quality needles that are gentle and easy to use. On the cover: You can see how the J-Type needle is especially gentle thanks to its extremely fine tip and shaft which allowed the needle to penetrate the balloon without bursting it!

In addition to conventional treatment methods in TCM (Traditional Chinese Medicine), you will find all our pioneering 3B LASER treatment products on page 24. The new and improved 3B LASER NEEDLE will make your work easier and with the improved cable guide, your cables will not become tangled and you'll have more time for your patients! The single-use applicators for 3B LASER NEEDLES are even more skin-friendly, more applicators than ever, and there are now 20 applicators more per pack.

Sit back, relax, and flick through the pages of our catalogue. We hope you enjoy reading about our products!

Kind Regards,



Enrique Dorn Schohaus
Acupuncture Product Manager

P.S. Looking for more? Please visit 3bscientific.com to view the complete product range, tutorials, videos, and much more.

More information about the gentle J-Type on page 6.

CONTENTS

SEIRIN® Product Range	4-11
Acupuncture Needles	
B-Type	4-5
J-Type	6
J-15	7
J-ProPak10	7
L-Type	10
M-Type	11
Gentle Press Needles	
NEW PYONEX	8
Spinex	11
Moxa & Accessories	12-16
MOXSAFE, Needle Moxa	16
Shōnishin	17
Cupping, Massage & Wellness	18-23
Cupping Instruments	18-19
Oils & Gels, TDP Lamps, Plasters	20
Dermarollers, Finger Rollers, Gua Sha	21
Massage Tools	22-23
3B LASER	24-29
3B LASER NEEDLE	24-26
3B LASER PEN	27
3B LASER SHOWER, Safety Glasses	28-29
Acupuncture Models	30
Wall Charts	31
Taping	32-33
Kinesiology Tape	32
Pre-Cut Tape	33
Massage Tables	34



SEIRIN® B-Type

Optimised single-use needles



B-Type - The best just got better!

For decades, SEIRIN Japan's high-quality B-Type acupuncture needles have been the world's leading single-use needles. The ergonomically optimised and precisely centred handle ensures safe and convenient use. The form, design and colour of the handle of the B-Type needle have been improved.

You will be able to get a better grip of the handle thanks to the new corrugated form. SEIRIN® is constantly working on improving their products even further with the goal of making the everyday job of the therapists easier - always keeping in mind the comfort of the patient. Your patients will thank you for it!

- Filed to diamond sharpness and polished three times
- No-slip grip with new corrugated form of handle
- Environmentally-friendly handle made of polypropylene plastic
- Electrolyte-cleaned
- Sterilised with EO gas
- Needle tip with extremely thin coating of medical silicone

► Use SEIRIN® needles for gentle, painless, and stress free acupuncture!

new



SEIRIN® B-Type/purple

Body acupuncture
Diameter x Needle length: 0.25 x 30 mm
Box of 100 needles
9935-1017768



SEIRIN® B-Type/purple

Body acupuncture
Diameter x Needle length: 0.25 x 40 mm
Box of 100 needles
9935-1017650



SEIRIN® B-Type/brown

Body acupuncture
Diameter x Needle length: 0.30 x 50 mm
Box of 100 needles
9935-1017653



SEIRIN® B-Type/black

Body acupuncture
Diameter x Needle length: 0.35 x 50 mm
Box of 100 needles
9935-1017654





Colour Coded handles for easy and fast size selection

SEIRIN® type B needles are available in 6 different gauges and lengths so the therapist can choose the right needle for different body sizes and areas.

The handles are colour coded to ensure ease of use. The colours of the new handle are now brighter and bolder, making them easier to see and faster to find.



Quick opening

The packaging can be opened easily and quickly with one hand. No more struggling to open packages! The different handle colour for each size of needle is indicated on the packaging allowing immediate recognition of the needle size. This is especially handy when you are working with different needle sizes.

Easy removal

The tried and tested easy-to-open blister pack can now be opened even more easily thanks to it being 5 mm longer. With this packaging, you can choose whether you want to use just one needle or all of the needles. Furthermore, SEIRIN® has improved the seal of the blister pack. There is no chance of impurity or contamination.



SEIRIN® B-Type/red

Cosmetic acupuncture

Diameter x Needle length: 0.16 x 15 mm

Box of 100 needles

9935-1017648



SEIRIN® B-Type/blue

Ear/body acupuncture

Diameter x Needle length: 0.20 x 15 mm

Box of 100 needles

9935-1017649



SEIRIN® B-Type/brown

Body acupuncture

Diameter x Needle length: 0.30 x 30 mm

Box of 100 needles

9935-1017652



Electrolyte-cleaned, sterilised with EO gas, and useable for at least 3 years.



For greatest precision, the needle has been incorporated exactly in the middle of the handle.



Rounded edges enable comfortable, painless positioning.

J-Type

J-Type – Virtually pain-free acupuncture
Acupuncture doesn't have to be painful. With the SEIRIN® J-Type needle, the treatment will be virtually painless for your patients. Due to new manufacturing methods, the needle tip and shaft are so fine that these SEIRIN® needles can be inserted with even less pressure and therefore especially gently by light tapping. They are made of special stainless steel and centred precisely in the handle, enabling particular insertion accuracy. The handle is very easy to grip due to its ridged surface, and is made of extremely lightweight polypropylene plastic. SEIRIN® type J needles are available in different gauges and lengths, which are colour coded.

Always on the safe side: SEIRIN® J-Type acupuncture needles are subjected to 8 stringent quality controls during manufacturing. They are electrolyte-cleaned, and sterilised with EO gas.

Box of 100 needles

Product Title / Colour	Application	Diameter x Needle length	Item No.
SEIRIN® J-Type / Dark Green	Cosmetic	0.12 x 30 mm	9935-1002412
SEIRIN® J-Type / Green	Body	0.14 x 30 mm	9935-1002414
SEIRIN® J-Type / Red	Body	0.16 x 30 mm	9935-1002416
SEIRIN® J-Type / Red	Body	0.16 x 40 mm	9935-1002417
SEIRIN® J-Type / Yellow	Body	0.18 x 30 mm	9935-1002418
SEIRIN® J-Type / Yellow	Body	0.18 x 40 mm	9935-1002419
SEIRIN® J-Type / Yellow	Body	0.18 x 50 mm	9935-1015548
SEIRIN® J-Type / Blue	Body	0.20 x 30 mm	9935-1002420
SEIRIN® J-Type / Blue	Body	0.20 x 40 mm	9935-1002421
SEIRIN® J-Type / Blue	Body	0.20 x 50 mm	9935-1002422
SEIRIN® J-Type / Blue	Body	0.20 x 60 mm	9935-1015549
SEIRIN® J-Type / Purple	Body	0.25 x 30 mm	9935-1002423
SEIRIN® J-Type / Purple	Body	0.25 x 40 mm	9935-1002424
SEIRIN® J-Type / Purple	Body	0.25 x 50 mm	9935-1002425
SEIRIN® J-Type / Purple	Body	0.25 x 60 mm	9935-1016664
SEIRIN® J-Type / Brown	Body	0.30 x 30 mm	9935-1002426
SEIRIN® J-Type / Brown	Body	0.30 x 40 mm	9935-1002427
SEIRIN® J-Type / Brown	Body	0.30 x 50 mm	9935-1002428
SEIRIN® J-Type / Brown	Body	0.30 x 60 mm	9935-1002429

J15-Type

J15-Type – For the most sensitive patients

J15-Type the thinnest and shortest J-type needle especially made for face/cosmetic acupuncture and are especially suitable for sensitive patients. J15-Type needles are 15 mm long and due to new manufacturing methods, the needle tip and shaft are so fine that these SEIRIN needles can be inserted with even less pressure and therefore especially gently by light tapping.



3B Scientific - Acupuncture Ears p. 31

Box of 100 needles

Product Title / Colour	Application	Diameter x Needle length	Item No.
SEIRIN® J-15 / Dark Blue	Cosmetic	0.10 x 15 mm	9935-1015547
SEIRIN® J-15 / Dark Green	Cosmetic	0.12 x 15 mm	9935-1002411
SEIRIN® J-15/ Green	Cosmetic	0.14 x 15 mm	9935-1002413
SEIRIN® J-15 / Red	Cosmetic	0.16 x 15 mm	9935-1002415
SEIRIN® J-15 / Yellow	Cosmetic	0.18 x 15 mm	9935-1017320

J-ProPak 10

J-ProPak10 – Functional and environmentally friendly

The J-ProPak10 by SEIRIN® is an intelligent development of the globally successful acupuncture needle. The particularly gentle new J-type has, innovative packaging design does away with individual guide tubes and therefore only needs half the amount of packaging material of the tried and tested J-type. That also means less waste. Each easy-to-open blister pack contains 10 acupuncture needles with a synthetic grip and only one guide tube.

So, this type of needle will be of interest to both fans and opponents of guide needles. During treatment, the blister pack turns into a functional needle dispenser, so that each needle remains sterile and yet is still easily accessible until it is used.

SEIRIN® J-ProPak10/blue - body Box with 100 needles


Diameter x Needle length: 0.20 x 30 mm

9935-1015552

SEIRIN® J-ProPak10/red - body Box with 100 needles

Diameter x Needle length: 0.16 x 30 mm

9935-1015551



Secure and long lasting application (up to 48 hours)


Breathable, skin friendly plaster that is almost invisible

NEW PYONEX

NEW PYONEX press needles guarantee gentle, superficial, and virtually imperceptible needle insertion. These are safe, easy-to-use, hygienic, and almost invisible at first sight.

Another advantage of this plaster is that parents can continue to treat their children gently at home.

The 0.20 mm thin needle is concealed within the adhesive surface of the plaster and can scarcely be felt by the patient.



Long lasting application

NEW PYONEX needles are embedded in a coloured plastic carrier pack, which does not look like medical packaging and is therefore especially popular with younger patients. The breathable plaster is stuck to the skin and can be fixed on precisely, safely and for long periods (up to 48 hours) and kept on even when washing and showering.




Easy to use at home

The needles are stimulated by pressing a finger down on the easily-to feel pressure point incorporated into the plaster. This ease of use helps make the patients be more involved in their own therapy.

Gentle acupuncture

NEW PYONEX needles can be used on a patient's body and, depending on the patient's anatomy, also on ear points. As the needles are available in short lengths (0.30-1.50 mm), NEW PYONEX gentle press needles are popular with children, patients with a low pain tolerance and patients with a needle phobia.





Box of 100 needles

Product Title / Colour	Application	Diameter x Needle length	Item No.
SEIRIN® NEW PYONEX / Orange	Ear/body/children	0.20 x 0.30 mm	9935-1002468
SEIRIN® NEW PYONEX / Yellow	Ear/body/children	0.20 x 0.60 mm	9935-1002471
SEIRIN® NEW PYONEX / Green	Body	0.20 x 0.90 mm	9935-1002465
SEIRIN® NEW PYONEX / Blue	Body	0.20 x 1.20 mm	9935-1002464
SEIRIN® NEW PYONEX / Pink	Body	0.20 x 1.50 mm	9935-1002469



Intensity of stimulation:
Low
Diameter x Needle length:
0.20 x 0.3 mm
9935-1002468



Intensity of stimulation:
Low - Medium
Diameter x Needle length:
0.20 x 0.6 mm
9935-1002471



Intensity of stimulation:
Medium
Diameter x Needle length:
0.20 x 0.9 mm
9935-1002465



Intensity of stimulation:
High
Diameter x Needle length:
0.20 x 1.2 mm
9935-1002464



Intensity of stimulation:
Highest
Diameter x Needle length:
0.20 x 1.5 mm
9935-1002469



L-Type

L-TYPE – The stainless steel needle with guide tube

Every SEIRIN® L-type needle is made entirely of stainless steel and is particularly suitable for electro-stimulation and moxibustion. They are supplied in a plastic guide tube sealed inside sterile blister packs. The guide tube helps make the L-type needle glide smoothly into the patient's skin and makes it easy for the therapist to make the treatment pain free.



Box of 100 needles

Product Title / Colour	Application	Diameter x Needle length	Item No.
SEIRIN® L-Type / Blue	Body	0.20 x 30 mm	9935-1002430
SEIRIN® L-Type / Blue	Body	0.20 x 40 mm	9935-1002431
SEIRIN® L-Type / Blue	Body	0.20 x 50 mm	9935-1015553
SEIRIN® L-Type / Blue	Body	0.20 x 60 mm	9935-1015554
SEIRIN® L-Type / Pink	Body	0.23 x 50 mm	9935-1015557
SEIRIN® L-Type / Purple	Body	0.25 x 30 mm	9935-1002432
SEIRIN® L-Type / Purple	Body	0.25 x 40 mm	9935-1002433
SEIRIN® L-Type / Purple	Body	0.25 x 50 mm	9935-1017795
SEIRIN® L-Type / Purple	Body	0.25 x 60 mm	9935-1002434
SEIRIN® L-Type / Brown	Body	0.30 x 30 mm	9935-1015558
SEIRIN® L-Type / Brown	Body	0.30 x 40 mm	9935-1015559
SEIRIN® L-Type / Brown	Body	0.30 x 60 mm	9935-1002435

new M-Type

The new SEIRIN® M-Type stainless steel needle

Longer handle for optimised grip

5 mm have been added to the handle, bringing it to a total of 25 mm. The longer handle made of light stainless steel provides the grip needed for precise treatments. The rounded edges of the handle increase comfort for the therapist. The needle itself has been manufactured to be placed exactly in the centre of the handle.

The ideal guide tube

The guide tube has been optimised and made thinner, so that the needle fits into it exactly. This guarantees precise handling of the needle. Additionally, the new material used for the guide tubes is softer, making it even easier to use.

Innovative and modern technology

Thanks to cutting-edge production technology, the needle and the plastic guide tube are being connected firmly and on the spot. This gives you the possibility to single-handedly and with just a light touch of the finger disconnect the needle from the guide tube.



Box of 100 needles

Product Title / Colour	Application	Diameter x Needle length	Item No.
SEIRIN® M-Type / Green	Ear/cosmetic	0.14 x 30 mm	9935-1018792
SEIRIN® M-Type / Green	Cosmetic / Body	0.14 x 40 mm	9935-1018793
SEIRIN® M-Type / Red	Ear/cosmetic	0.16 x 30 mm	9935-1018794
SEIRIN® M-Type / Red	Body	0.16 x 40 mm	9935-1018795
SEIRIN® M-Type / Red	Body	0.16 x 50 mm	9935-1018796
SEIRIN® M-Type / Yellow	Ear/cosmetic	0.18 x 30 mm	9935-1018797
SEIRIN® M-Type / Yellow	Body	0.18 x 40 mm	9935-1018798
SEIRIN® M-Type / Yellow	Body	0.18 x 50 mm	9935-1018799
SEIRIN® M-Type / Blue	Ear/cosmetic	0.20 x 30 mm	9935-1018800
SEIRIN® M-Type / Blue	Body	0.20 x 40 mm	9935-1018801
SEIRIN® M-Type / Blue	Body	0.20 x 50 mm	9935-1018802
SEIRIN® M-Type / Purple	Body	0.25 x 30 mm	9935-1018803
SEIRIN® M-Type / Purple	Body	0.25 x 40 mm	9935-1018804
SEIRIN® M-Type / Purple	Body	0.25 x 50 mm	9935-1018805

SPINEX

SPINEX – Intradermal acupuncture needle

This very special needle was invented by Kohei Akabane in Japan. It is applied in four steps:

1. A tender point or an induration is found in the skin
2. The skin is stretched
3. The needle is inserted 2–3 mm into the dermis
4. The needle is covered with an adhesive tape and left in place for a timed period

Whilst remaining in the dermis, the acupuncture needle can provide powerful pain relief. Packaging: One sterile needle per blister in strips of 10. Each box contains 100 intradermal needles. The Spinex needle can be selected as required to achieve various stimulus strengths.



Box of 100 needles


Description	Application	Diameter x Needle length	Item No.
SEIRIN® SPINEX	Body	0.12 x 3 mm	9935-1002461
SEIRIN® SPINEX	Body	0.12 x 4 mm	9935-1002462
SEIRIN® SPINEX	Body	0.14 x 6 mm	9935-1002463

Moxibustion


Moxibustion is, just as with acupuncture, a tried and tested healing method in Traditional Chinese Medicine (TCM) and has been used for millennia.

Dampness and cold are removed from the meridians through the strength of the yang. Moxibustion is especially beneficial during the cold season.

The pleasant heat restimulates the flow of Qi, or life energy. Direct moxibustion is mainly used in China, while in Western countries indirect moxibustion is much preferred.



A



B

A. Gold Mountain Moxa

Made from carefully selected superior quality mugwort gathered only during mid-winter when it is at its peak quality from Ibuki Mountain. Super pure, extra fine texture make this the highest quality moxa available. Excellent for direct moxibustion. Pure moxa is less likely to burn the skin because of its high concentration of essential oil. Gold mountain has a volume limit due to its limited gathering and production period.

100 g 9935-1015604

B. Wakakusa Moxa

Wakakusa is the choice moxa for use in moxa needle heating (fine needle). Excellent adhesiveness with no stems. Its heat and easy-lighting make it ideal for applying to acupuncture needles.

30 g 9935-1015601



Arirang – Self Adhesive Moxa
 Burning time: 2 to 4 minutes; made in Korea.
 200 per box.
 9935-1002442



Red Mini Moxa
 Self adhesive moxa; intense heat 60°C;
 Burning time: 8 minutes; made in Korea.
 180 per box.
 9935-1002448



Green Mini-Moxa
 Self adhesive moxa; intense heat 60°C;
 Burning time: 8 minutes; made in Korea.
 180 per box.
 9935-1002449



Nien-Ying Moxa
 Contains mainly moxa. Mixed with a variety of herbs (Rhizoma Alpiniae Officinarum, Ramulus Cinnamomi, Radix Angelicae Pubescentis, Radix Angelicae Dahuricae).
 9935-1002452



Tai-Yi Moxa
 Contains mainly moxa. Mixed with a variety of herbs (Rhizoma Et Radix Notopterygii, Radix Aucklandiae, Radix Angelicae Dahuricae, Radix Salviae).
 9935-1015593



Pure Moxa Roll
 Classical Chinese pure moxa rolls. For indirect moxa treatment. Long and thick rolls with pure moxa contents.
 9935-1002451



A. Smokeless Moxa Perforated Rolls
 Charcoaled moxa cigar with a central hole for easy lighting.
 9935-1015598

B. Hoist Smokeless Mox
 Moxa rolls with less smoke. Very condensed charcoaled moxa offers slow burning and gentle heat. Contains a variety of herbs, which may generate an odour. 5pcs/pack.
 9935-1015600



Smokeless, Self Adhesive Mini-Moxa
 Burning time: 7 to 10 minutes; made in Korea.
 180 per box.
 9935-1002445



Moxa Sticks – Moxa Cigars
 Burning time: 11–12 minutes; made in Korea. 20 per box.
 9935-1002443

High quality moxa, smokeless and odor-free



A. Smokeless Gu Gwan Moxa

Excellent moxa quality for needle warming. Smokeless and odour-free. Charcoaled moxa rolled into a tube with hole in the middle. It is to be used with the metal cap (1015581 Needle Cap).

Insert the metal cap into the tube, place the metal cap on top of the acupuncture needle. It can also be used with moxa devices (1015589) and moxa burning bowls (1015585 and 1015586).

10pcs of metal caps
9935-1015582

B. Gu Gwan Moxa – 25 mm

Compressed moxa formed into small tube with a hollow centre in the middle. Very easy to use, whether holding it or with moxa holders 1015585 or 1015586 (page 13).

7mm diameter x 20mm
9935-1015584



Metal Cap for Gu Gwan Moxa

Insert the cap into the needle moxa (1015583 or 1015584), place it on the top of the needle handle. Suitable for SEIRIN L-Type needles (page 10). Not suitable for Chinese style needles as they have a loop on the handle. Reusable.
9935-1015581



Gu Gwan Moxa - Long 105 mm

Compressed moxa formed into small tube with a hollow centre. Very easy to use, either holding it or with moxa holders 1015585 or 1015586.

7 mm diameter x 105 mm
9935-1015583



Holder for Gu Gwan Moxa

Ideal for indirect moxa treatment. Place the device on the treatment area with moxa sticks burning. Holes on the device allow acupuncture needles through.
9935-1015585



Moxa Device

Screw type device for moxa treatment. It has mesh in the middle of the device to hold the burning moxa, the heat passes up and down through the mesh protected top and bottom. The screw in the middle can adjust the height of the device controlling level of heat.

Secure one device, or more if necessary, with the belt and then attach around the waist or abdomen. Also recommended to use with moxa pads for large areas such as shoulders and back.

9935-1015589



Charcoal Moxa Cone

A cone shaped moxa generates very gentle and longlasting heat for moxibustion treatment. The unique oven baking technique produces moxa charcoal which eliminates smoke, tar, and odour associated with traditional moxa. The cone is coated with a special compound comprised of Bio-clay, Jade, Germanium and Bio-ceramic for enhanced effect.

It requires approximately 5 minutes until effective temperature is reached. Each cone lasts for about 20 to 25 minutes, the residual heat is still powerful when complete. Highly recommended to use this product with DONGBANG Moxa Device.

9935-1015588

Ondan Smokeless Needle Moxa

120pcs with 1 scoop and 6 basket holders.

Japanese odour-free smokeless moxa. 12mm long cylinder shape with a 3mm hole in the middle to go through the needle handle. Cylinder diameter is about 12mm.


9935-1015594



Basket for Smokeless Needle Moxa

Place the basket on top of the needle handle, place the smokeless moxa (1015594) on top of the basket.

9935-1015595



Equipment for safe moxibustion treatments!



Tiger Warmer

Tiger warmer device for powerful heat on larger areas. Excellent heat transmitting and protection.

9935-1002446



Lion Warmer Device

Burn the moxa inside and warm acupoints along the meridians. Excellent heat transmittance and protection. To be used with Long Gu Gwan moxa (1015583).

9935-1015587



Smokeless Moxa Perforated Rolls

Charcoaled moxa cigar with a central hole for easy lighting. 5 mm diameter

9935-1015599



Tiger Warmer Moxa

Compressed moxa stick for Tiger Warmer device. Thin but solid. Size: 5mm diameter x 80mm Length. 20 sticks per box. Made in Korea

9935-1002447



Moxa Spoon

Must-have for indirect moxibustion: the Moxa Spoon. Stainless steel spoon with a split in the middle. Safe and easy to remove moxa from acupuncture needle.

9935-1015592



Moxa Roll Extinguisher with Plate

Moxa extinguisher mounted in the middle of the plate. Safe and stable to extinguish, hold and collect the moxa rolls.

9935-1015590

MOXSAFE

Safe-smokeless-odourless needle

MOXSAFE, which is produced in Japan, is the newest development in moxibustion. Its name says it all – this type of moxibustion is safe and produces very little smoke.

MOXSAFE consists of a stainless steel ring-shaped moxa carrier and a mugwort carbon ring.

MOXSAFE is placed on the metal needle and easily lit by a lighter. The heat produced is transferred directly to the needle and thus penetrates the adjacent regions of skin. The mugwort charcoal burns for exactly 6 minutes. The solid ash remains safely on the moxa carrier. After a cooling-down phase, the MOXSAFE is removed and disposed of with the metal needle.

A. SEIRIN® Moxsafe 20

Pack with 20 moxa and 20 moxa carriers
9935-1002457

B. SEIRIN® Moxsafe 100

Refill pack with 100 moxa and 6 moxa carriers
9935-1009707



The perfect combination!

When MOXSAFE is used with L-Type needles (min. 30 mm in length, see page 10), it provides a safe and precise heat treatment for therapists and patients.



Needle Moxa

Burning time: 12–15 minutes; made in China.
200 per box
9935-1002444



OS Moxa Balm

Ingredients: moxa, herbs, etc. Multi-functional: can be used for direct moxa, infrared heating treatment, and massage.
9935-1015596

Yoneyama Needle for Children

A good general purpose tool for treating children. Used for rubbing or scratching the skin for cutaneous treatment along meridian pathways or with special techniques. All of the surfaces and corners may be used to provide a different stimulation. 36 (L) x 21 (W) mm
9935-1015647



Tsun Meter Point Maker


Interprets 1 Tsun into different distances (0.5, 1.5, 2 Tsun) to easily locate acupuncture points. One Tsun is the distance between the creases of the middle finger.
9935-1015652



Hygienic and safe!

Sharps Box

The new disposal box provides an absolutely safe method of properly collecting and disposing of used acupuncture needles. Of course, all boxes fulfil the most recent international standards for disposal boxes. Box capacity is between 3,000 and 5,000 needles, depending on the type and size of the needles.
9935-1004882



A. Kakibari Scratching Needle for Children

This triangle shaped needle can be used more for individual point stimulation along the meridian. Either end may be used for scratching or stimulating. Used for the treatment of children. 36 (L) x 6 (W) mm.
9935-1015648

B. Choki Shōnishin Needle for Children


Every angle of this instrument can be used for varying intensities of cutaneous treatment. 80 (L) x 30 (W) mm
9935-1015649



C. Ear Point Finder - simple
9935-1002480


D. Ear Point Finder - double-sided
9935-1002481

A must-have for every practice!



Meridian Notepad

Our popular, high-quality A4 wall charts; diagram of the key body and ear points; multicoloured; improved patient documentation; 50 pages.
9935-1002440



Tsumo-Shin Spring Modulator with Magnet

The most popular, adjustable pressure, spring-activated cutaneous probe for treating or locating trigger points, acupoints, or ear points.

5" long modulator which includes 3" x 0.4" diameter grooved handle and 1.8" spring loaded probe with rounded tip (does not puncture skin). The other end features a powerful 2,000 gauss magnetic ball probe (4 mm diameter) which may also be used for acu-pressure and acu-massage.

The polarity may be easily changed from bio-north to bio-south pole by removing the ball probe and reversing the magnet's position. Includes five interchangeable, different strength springs:

- 100 g = Green
 - 150 g = Yellow
 - 200 g = White
 - 250 g = Purple
 - 300 g = Red
- 9935-1015650



image similar



MERIDIUS Cupping Set
Ideal for in-home visits!

The newly introduced MERIDIUS cupping set include: 17 cups, 1 hand pistol pump and 1 connecting tube, conveniently stored in a navy blue case, ideal for travelling. The blue silicone valve can be removed for sterilising the jar. Versatile thanks to the various jar forms. Thanks to the thicker material and the solid workmanship, the jars are very robust and long-lasting.

9935-1015606



Plastic Cupping Set (17 cups)

Elegant lightweight case, easy to carry and store. Cups are non-breakable and there is a valve on the top. Place the hand pump on top of the cup, apply pressure. To release pressure press the valve.

- 17 cups in total
- 2" diameter (13pcs)
- 1.8" (1pc)
- 1.6" (1pc)
- 1.3" (1pc), and 1.0" (1pc)
- 1 suction pump
- 1 extension tube, 26" length

9935-1013327

Fire cupping glass set (not shown)

Cupping glasses are made from heat-resistant borosilicate glass and can be easily cleaned and disinfected (dishwasher, steriliser etc.). These high-quality cupping glasses are made by hand (mouth-blown and hand-shaped).

Delivery contents: 6 cupping glasses, approx. 40 mm

40 mm 9935-1018093

50 mm 9935-1018094

60 mm 9935-1018095

70 mm 9935-1018100

▶ Made by hand for highest quality!



DB Disposable Plastic Cups

Disposable plastic cups for single use.
Sterile. Compatible with 1015607 hand pistol pump.
5 sizes available. 5pcs/pack.

Size	Quantity	Item No.
5.6 cm	5	9935-1015609
5.1 cm	5	9935-1015610
4.6 cm	5	9935-1015611
4.0 cm	5	9935-1015612
3.3 cm	5	9935-1015613



Plastic Cupping Set (10 cups)

- 10 cups in total
 - 2" diameter (6pcs), 1.8" diameter (1pc), 1.6" diameter (1pc), 1.3" diameter (1pc) and 1.0" diameter
 - No extension tube included in this set, can be purchased separately
- 9935-1015605**



Massage Glass Cupping Set with Suction Bulb

Medical massage-glass set with suction ball for massaging and fire-free cupping. The massage-glass set is made of glass so you can clean and disinfect it easily (dishwasher safe). They are handmade in the town of Lauscha, Germany, to a very high standard.

Size	Quantity	Item No.
40-60 mm	3	9935-1017757
25-45 mm	5	9935-1017754
20-60 mm	5	9935-1017755
30-70 mm	5	9935-1017756
20 mm	6	9935-1018620
40 mm	6	9935-1018098
50 mm	6	9935-1018090
60 mm	6	9935-1018091
70 mm	6	9935-1018092



Rubber Cupping Jars

- Set of 4
 - 70 mm, 60 mm, 45 mm and 30 mm diameter
- 9935-1015608**



Hand Pistol Pump

Replacement hand pump for applying pressure with cupping therapy. Compatible with disposable cups.
9935-1015607

Lian Ru Emulsion / Gels

Lian Ru relax (warming)

Lian Ru relax (warming) promotes the flow of Qi and blood, and relaxes and warms the muscles and tendons. It is ideal for tension and muscle pain with a feeling of coldness, especially if the problems worsen with damp, cold weather.

75 ml 9935-1015653

Lian Ru (cooling)

Lian Ru (cooling) promotes the flow of Qi and blood, and can be used in sports as an ice spray or ice pack.

75 ml 9935-1015656

Lian Ru (relax)

When the body then needs to regenerate after exercising, Lian Ru relax promotes the flow of Qi and blood, frees up the pathways and relaxes muscles and tendons. Lian Ru relax slightly warms the muscles and is perfect for when there are no signs of distinct heat or cold sensations.

75 ml 9935-1015657

Lian Shui Fluid Spray

Lian Shui Fluid Spray - relaxes, warms and regenerates the muscles and stimulates Qi and blood. Lian Shui is quickly absorbed into the skin, bestowing a feeling of freshness. Activates Qi and blood, removes stagnation and stasis, has a slight warming effect and opens up the pathways. artificial colours and preservatives.

100 ml 9935-1015659

Lian You Massage Oil

Ideal for massaging, cupping and for Gua Sha treatment. Relaxing massage oil with an extract from soothing Chinese herbs and quality plant oils. A fresh fragrance composition made from 100% natural essential oils of sandalwood, ginger, incense and lemongrass stimulates the senses.

50 ml 9935-1015660

500 ml 9935-1015661



▶ Long lasting warming effect with 100% natural ingredients!



HerbaChaud® TDP Heating Lamp, Modell 3000 with Timer

The patented HerbaChaud® TDP-electro magnetic thermal therapy lamp is an ideal enhancement of your treatment and therapy options. It was developed in China by a group of scientists and physician under the leadership of Dr. Gou. In China it ran through hundreds of tests in research institutes, clinics and universities. Unlike conventional therapeutic infrared and heat lamps, the HerbaChaud® TDP-Lamp has a round sheet, coated with **33 natural minerals**, which has a soothing and healing effect on the human body.

CE/TÜV approved, registered with Swissmedic, digital timer, replacement sheet for additional 1,000 operating hours available.

9935-1018806



HerbaChaud® Heating Plasters

The core of the HerbaChaud® Heating plaster consists of a traditional mix of Chinese medicinal herbs, enveloped in a sachet made of orange tree bark. On the side facing the body, it is covered with a black layer composed of a carbon-iron mixture. This layer is covered with minerals that are vital for the human organism. The heat develops as iron oxidises upon exposure to air and emanates in a controlled way, providing constant and pleasant warmth for a period of more than 12 hours.

Advantages of the HerbaChaud® Heating plaster:

- Constant heating effect
- Kind to the skin, no itching or burning
- Revitalising and relaxing effect
- HerbaChaud® heating plaster is made of 100% natural and plant agents
- Very easy to use
- Active substances contained in the HerbaChaud® Medicinal herb mixture
- Mugwort, Ginger, Chinese cinnamon, Angelica sinensis, Notoptergium root
- Costus root oil, Carbon-iron mix, Mineral mix

Pack of 2 9935-1005927

Natural Heating Patch

The patch's main effect is achieved by the generation of heat, which is produced when the iron inside the patch oxidises on contact with the air, and is controlled to ensure a consistent and pleasant warming effect for 16 to 20 hours.

- 100% natural
- Kind to the skin
- No side effects, even when used for prolonged periods of treatment
- Soothes aches and pains
- Stops inflammation
- Relaxes and revitalises

9935-1017314



Hydrogel Cooling Plaster

Cooling therapy for sports injuries and muscle pain is well established and widely accepted. Prompt cooling has been shown to effectively reduce pain. The hydrogel plaster works like ice by cooling the applied area. Cooling therapy can also be used to soothe tense or stiff muscle. Hydrogel cooling plaster is convenient and easy to apply, can be placed directly onto the skin and used under clothes. For best result apply the patch to the required area. The cooling effect from the patch can be felt quickly. Generally, one patch can last up to 8 hours. Repeated applications may be required. Suitable for daily use.

9935-1017315

Derma-Q-Roller

The DERMA-Q Microneedle Roller is used to reduce wrinkles, scars, stretch marks, cellulite, and hair loss. Derma-Q rollers are ideal for perioral wrinkling, stretch marks, fine lines and acne scarring. Clinical studies have shown microneedling therapy is effective.

The creation of micro-channels within the dermal layers has shown to stimulate collagen, increasing production levels to thicken skin, erase wrinkles, and smooth scars.


Application areas are face, neck, and body. Avoid eye area. Available in 5 sizes: 0.2mm, 0.5mm are recommended for beginners, 1.0mm is recommended for experienced users, 1.5mm, 2.0mm are recommended for experienced or professional users.

Derma-Q Classic, 192 Needles.
Derma-Q Premium, 260 Needles, gold-plated.

Derma-Q Premium


Derma-Q Classic

Description	Length of needles	Item No.
Derma-Q Classic	0.2 mm	9935-1015624
Derma-Q Classic	0.5 mm	9935-1015625
Derma-Q Classic	1.0 mm	9935-1015626
Derma-Q Classic	1.5 mm	9935-1015627
Derma-Q Classic	2.0 mm	9935-1015628
Derma-Q Premium	0.2 mm	9935-1015629
Derma-Q Premium	0.5 mm	9935-1015630
Derma-Q Premium	1.0 mm	9935-1015631
Derma-Q Premium	1.5 mm	9935-1015632
Derma-Q Premium	2.0 mm	9935-1015633



Finger Roller

5 rollers per pack.
 Gold-coloured (Yang) 9935-1002484
 Silver-coloured (Yang) 9935-1002485



Spring loaded dermal roller with probe

Multi-functional, dermal roller with spring active probe. Silver-plated for sedating purposes. Rolling heads with pyramid shaped surface. The shaft can also be used to massage hand and foot.
 9935-1002482

Smokeless Incense

Sandlewood scented, smokeless incense, for lighting moxa. Length: 170mm. 200pcs / box
 9935-1015622



Gua Sha Jade Stones

Jade stones for Gua Sha massage.
 Classic 9935-1015615
 Heart-shaped 9935-1015614



Premium Gua Sha Set

- One roller
 - One hammer shape
 - One boomerang shape
 - One jade classic shape
 - One horn shape
 - Two comb shape
 - One bottle of treatment oil
- 9935-1015616



Ceramic Wonder

► Can be used warm and cold!

Now you can experience the **next best thing to a real massage** in your own home, at work or just about anywhere. For more therapeutic effect, the Touch Down personal massage tool can be used heated or chilled. Simply warm up the tool under hot water or chill in a refrigerator or in ice. Elegantly designed, the therapeutic Ceramic Wonder Touch Down massage tool can be used to ease muscle pain and stress, providing relief to practically everyone—from children to seniors and from the competitive athlete to the busy executive. This fantastic **eco friendly** device is designed to provide you with a **muscle-relaxing, stress-reducing, stimulating massage experience**.

The Ceramic Wonder Touch Down can be used in many ways to apply deep pressure to muscular trigger points and tension or with semi-circle moves over the rest of your body. It's great with oil or body creams.

Pamper your Pets!

Even your pets will love a massage with the Ceramic Wonder Touch Down!

Ceramic Wonder Touch Down Massage Tool Original

9935-1015666

Ceramic Wonder Touch Down Massage Tool Mini

9935-1015667



Ceramic Wonder Professional Tool Set, 2 pieces

The Ceramic Wonder Professional Tools Set is a **hand-crafted massage** tool made from natural and environmentally safe clay and glaze. This set includes one 9" roller and one 14" roller. Rollers are **lightweight and durable**. They aid in reducing therapist fatigue. Heated, the rollers serve as an extension of the therapist's hands to warm and loosen tight muscles and balance energy points in the body. Chilled, the rollers are a powerful decongestant as it pushes fluid and waste from overworked, tense or inflamed tissues. Use the end of rollers for pressure point therapy. Includes cloth storage bags with drawstring closure. **9935-1015665**



Ceramic Wonder Spine Fork Tool

This tool is ideal for working along the erector spinae (extensor spinae) muscles. The design of the Spine Fork allows the user to work on very specific areas along the spinal column. The tool can also be very useful when working on the entire back by gripping the tool on both ends and using the centre or entire body of the tool to perform a variety of strokes. Lightweight and durable, this can be used heated or chilled.
9935-1015668



Ceramic Wonder Thumbs Up Tool

This tool is designed to aid the therapist's thumb. Lightweight and durable, this can be used heated or chilled.
9935-1015669



Ceramic Wonder Reflex Globes, Set of 2

- Relieve tension
- Promote relaxation
- Restore energy to the face and body
- Increase circulation to promote oxygenation of the skin
- Stimulate collagen fibres to improve elasticity
- Tighten and lift the facial muscles

Made of eco friendly clay and glaze, these light weight and durable globes can be used heated or chilled.

Simply warm up the globes under hot water or chill in a refrigerator or ice. Includes storage bag with drawstrings.
9935-1015670



A. Massage Cube

Traditional Asian massage cubes made of polished hardwood can be used to specifically treat trigger points with little effort, and work right into the deep muscles.
9935-1015840

B. Massage Rod

The classic in trigger point therapy. By placing pressure on specific areas, blockages and tensions can be loosened up and removed. The massage rods can be used for a particularly deep massage of the muscles.
9935-1015843

C. Massage Birdy

The massage birdy, thanks to its shape, fits snugly in the hand. It allows the targeted treatment of trigger points deep into the muscles by applying specific effort.
9935-1015842

D. Massage Cross

Like the massage Birdy, the cross form fits especially well in the hand and reduces the effort required by the therapist.
9935-1015841





3B LASER NEEDLE

compact | economic | simple to use

OFFER YOUR PATIENTS A PAIN-FREE ACUPUNCTURE TREATMENT

- Treat up to 12 acupuncture points
- Space-saving thanks to compact structure
- Simple to use
- Low life cycle costs

The innovative 3B LASER Therapy is classified as low-level laser therapy, and is a **painless and noninvasive treatment procedure**. Light particles (photons) with high excitation density are conducted through fine optical fibres. The fibres do not puncture the skin, but are placed on the therapy points, radiating into the tissue which requires treatment. Laser needles

can be used simultaneously on any point combinations on the body, in the head region and on the ear. Thus, corresponding methodically to the traditions of classical acupuncture. They are equivalent to acupuncture treatments in their cerebral and peripheral stimulation effects. Additionally, the laser beams stimulate complex healing processes in tissues.

You control and monitor the 3B LASER NEEDLE via the **mobile touch screen** (Control Touch Pad), which is connected with the laser unit by a radio interface. As a result of this novel technology, you have more flexibility in your practice routine as a therapist. Additionally, you achieve optimum communication with your patients, who can follow the treatment at all times and thus be more actively involved in the elucidation and therapy. This creates additional trust and increases the therapy success. The **daily work is greatly simplified** so that both a bedside and remote observation of the therapy course is possible. The input and control of the treatment parameters is self-explanatory and can be easily delegated to qualified staff members.

In addition to the preset programs and frequencies (Multi, Alpha, Nogier, Bahr, Reininger), **individual patient programs** can be easily entered, stored under the patient's name and retrieved. In this manner, you can incidentally prepare your own database. Power (mW), frequencies (Hz), treatment time (min) and number of laser diodes can be varied.

Product Specifications:

- Classification Class IIa according to EU Guideline RL 93 / 42 EC
- Laser class 3R according to EN 60825-1:2008
- Effective output in CW operation 12x50 mW (1004916); adjustable in steps of 1 mW
- Number of laser needles: 12
- Wavelength 660 nm; additional wavelengths on request
- Modulation frequencies 1-10,000 Hz
- Conformity CE 0366

For more information and additional products please look at the following pages!



3B LASER NEEDLE

With the new 3B LASER you use up to **12 lasers simultaneously**. An instrument version with red laser light (660 nm) is provided as a proven standard. The instruments can be individualised according to the focus of your therapy, so that in addition to red laser light (660 nm), infrared light (808 nm) and blue laser light (405 nm) can be used and combined.

All advantages at a glance:

- 3B LASER unit with 12 laser diodes
- Mobile Touch Screen (Control Touch Pad)
- Self-explanatory and logical operation
- Low space requirements and high level of mobility
- Perfect Ergonomics
- Reliable and sophisticated technology
- Painless and side-effect-free treatment
- Broad therapy spectrum
- High patient acceptance
- Low capital expenditure
- Pays for itself after only a few treatments
- Comprehensive, fast service
- Made in Germany: guaranteed highest product quality and safety
- Multilingual Product Manual (EN, FR, ES, IT, PT, DE)

12 x 660 nm (red laser light), 50 mW, in total 600 mW
9935-1004916

10 x 785 nm (infrared laser light) and 2 x 660 nm (red laser light), 50 mW, in total 600 mW
9935-1008814


► For more information please visit 3bscientific.com



Optimised functionality

Via the colour display of the touch screen, you are always informed of the important functions of the instrument, such as treatment time and dosage.

In addition to the preset programs and frequencies (Multi, Alpha, Nogier, Bahr, Reiningger), individual patient programs can be easily entered, stored under the patient's name and retrieved. In this manner, you can incidentally prepare your own database. Power (mW), frequencies (Hz), treatment time (min) and number of laser diodes can be varied.



Flexibility when you need it

An extremely flexible treatment arm and intelligent cable routing ensure very good handling and optimum treatment. Adhesive applicators, which simplify affixing the waveguide and save the therapist valuable time, are provided as standard.



3B LASER NEEDLE Fibre-Optic Cable Holder Extension
9935-1019258



3B LASER NEEDLE Head Adaptor for Ear and Scalp Acupuncture
The head adaptor enables the therapist to reach the acupuncture points in the head area safely and reliably, in particular for ear acupuncture. The laser needles can be placed with precision thanks to the head adaptor
9935-1012412



3B LASER NEEDLE Dermatology Spot
The flexible 3B Dermatology Spot is a useful supplement to the 3B LASER NEEDLE (1004916 and 1008814). With the flexible Dermatology Spot it is possible to treat larger skin areas individually.

Applications: Dermatology, Orthopaedics, Sports Medicine, Physiotherapy.
9935-1004905



new

3B LASER NEEDLE Single-Use Applicators
The disposable applicators for 3B LASER NEEDLE are made of high quality plastic. The acrylic adhesive is hypoallergenic and leaves no residue on the skin. The disposable applicators are hygienic and have been developed specifically for the 3B LASER NEEDLE. One pack contains 120 applicators.
9935-1018635



Side Table for 3B LASER NEEDLE
Practical side table for 3B LASER NEEDLE.
Dimensions: 60 cm x 40 cm x 90 cm (W x D x H).
9935-1009656



Wall Mount for 3B LASER NEEDLE
The 3B LASER NEEDLE wall mount is the perfect solution for small rooms
9935-1009657



Cart for 3B LASER NEEDLE
The mobile single-column cart for 3B LASER NEEDLES features a slim design and despite the small footprint, this equipment card can take a load of up to 65 kg.
9935-1017796





Table mount for 3B LASER NEEDLE
The universal table bracket is suitable for all tables.
9935-1009658

3B LASER NEEDLE Optical Fibre for Infrared Light (not shown)
Flexible optical fibre for infrared laser light.
Replacement part for 3B LASER 1017714 & 1008814, white; 2.5 m.
9935-1013292

3B LASER NEEDLE Optical Fibre for red light (not shown)
Flexible optical fibre for 660 nm red laser light.
Replacement part for 3B LASER 1004916.
9935-1008813



3B LASER PEN

ergonomic and powerful

With its 200 mW, the 3B LASER PEN series is 50 mW more powerful than its predecessor, which means treatment times are shortened, and the number of indications is even higher. You can set the treatment time or the dose, and the 3B LASER PEN 200 mW will calculate the rest automatically. This makes your everyday work easier.

In addition to the pre-set frequencies (Nogier, Bahr, Reininger, Chakra), you can programme frequencies yourself quickly from 0.1 to 10,000 Hz, and then store these settings (20 storage locations). A practical charge station rounds off the features of the 3B LASER PEN 200 mW and makes it one of the forerunners of its class.

Product features:


- Display for left and right-handed operators
- Information on therapy time, laser performance dose, frequency and charge status
- Memory function for last values
- Integrated point finder
- Adjustable laser performance
- Wavelength 808 nm, infrared laser light
- CW laser beam
- Nogier frequencies
- Bahr and Chakra frequencies
- Reininger frequencies
- Programmable frequencies
- The treatment time or dose can be set

All lasers come complete with safety goggles for both the therapist and patient, an acrylic stand and a charger in an aluminium case.

3B LASER PEN 200 mW, 808 nm
9935-1013278

3B LASER PEN 100 mW, 785 nm
9935-1015662

3B LASER PEN 100 mW, 660 nm
9935-1016649



A. Dental Wave Guide

Dental wave guide for 3B LASER PENS 100 mW and 200 mW
9935-1004889

B. Straight Wave Guide

Dental wave guide for 3B LASER PENS 100 mW and 200 mW
9935-1004890



Other Accessories:

Battery Cover for 3B LASER PEN 100, S100 and 200 mW

- 1013278 - 3B LASER PEN 200 mW, 808nm, Infrared
- 1016649 - 3B LASER PEN 100 mW, 660nm, red
- 1015662 - 3B LASER PEN 100 mW, 785nm, Infrared

9935-1017582

Set of 4 rechargeable batteries for all 3B LASERS
Replacement batteries for the 3B LASER PEN 100 and 200 mW and for the Control Touchpad of the 3B LASER NEEDLE
9935-1004898

Laser Warning Sign (not shown)
9935-1004899



LASER SHOWER POWER TWIN 21

Innovation for your practice. The new Power Twin 21 laser shower – holder of the reddot design award – has been specially designed for use on large areas and now offers **new treatment options**. For the first time, the laser shower and the laser comb have been brought together in a single device:

- More laser power thanks to 21 laser diodes
- Large beam area of 55 cm²
- Individually adjustable settings via the colour display
- Switch to a laser comb at the click of a button
- Very hairy skin can be treated
- Continuous ray and multi-frequency mode: Nogier, Bahr, Reiningger frequencies have been pre-set
- Laser power can be set at between 5 mW and 50mW

Technical information:

- CE 0366, Laser class: 3B
 - Total diode power: 1,050 mW (21x50 mW)
 - Wavelength: 785 nm (infrared)
 - Four red light diodes limit the treatment area
- Other power levels and wavelengths on request.

Includes acrylic storage bag and aluminium box.
9935-1004911

new

LASER SHOWER POWER TWIN 21; 21 x 100 mW

Higher performance: total diode output: 2100 mW (21x100 mW). Benefit from the doubled output and optimise your treatment successes.

9935-1018742



Laser Comb Attachment

Replacement for 1004911 - LASER SHOWER Power Twin 21.

9935-1004914



Patient Laser Safety Goggles

The laser safety goggles for the patient. To be used for both both red and infrared laser light.

9935-1004888



new

Safety glasses for infrared laser (785-808nm)

Speical features:

Application: Excimer, UV, harmonics of Nd: YAG

(266 nm), Ti: Sa, diode, disk, Nd: YAG, fibre

Filter material: Special synthetic material

Filter colour: Yellow greenish

Filter thickness: Approx. 2.2 mm

VL: Approx. 60%

Frame: TERMINATOR

Colour: Yellow / black

Certified protection levels according to EN 207:

190-315 D LB7 + IR LB4

760-790 DIR LB2

790-805 DIR LB3

805-840 DIR LB4

840-860 DIR LB5

860-1030 D LB5 + IR LB6

1030-1065 D LB5 + IR LB7

1065-1080 DIR LB5

1080-1100 DIR LB3

PF S CE

9935-1018738



Eyeball Glasses for all Laser Lights

Silver, high-grade steel; non-reflecting surface; protection levels from 315 nm–1400 nm

9935-1004912

3B LASER SHOWER

The 3B LASER SHOWER was specifically designed for extensive applications of laser therapy. The concurrent operation of 14 laser diodes allows an area of 35 cm² to be treated. The laser diodes are arranged in such a way as to enable the best possible therapeutic effects. Four intensively radiating red light diodes serve as guide beams and limit the area of treatment. A special filter glass protects the laser diodes and can be easily and hygienically cleaned. The 3B LASER SHOWER is a lightweight, portable and compact instrument. Simple and easy to use.


Applications:

- Dermatology
- Orthopaedics
- Sports medicine
- Physiotherapy

Technical data:

- Total diode performance: 700 mW (14 x 50 mW)
- Wavelength: 785 nm (infrared)
- CW (Continuous Wave) and multiple frequency operation

9935-1004901



new

Safety glasses for infrared laser (785-808nm)

Material: Metal

Colour: Titanium

Special features:

- Metal frame in noble titanium-look
- Wrap-around lens provides perfect side impact protection without constricting the field of vision
- Nose softpads (best comfort)
- Thin slots at the upper edge prevent fogging
- Light weight

9935-1018739



new

Safety glasses for red laser (660nm)

Speical features:

- sporty design
- Best fit via "wrap-around" lens
- Additional side protection
- Comfortable soft nose-pad
- Featherweight (only 26 g)
- Soft headstrap included in delivery

9935-1018741



new

Safety glasses for red and infrared laser (660 & 785-808nm)

Special features:

- Sporty design
- Best fit via "wrap-around" lens
- Additional side protection
- Comfortable soft nose-pad
- Featherweight (only 26 g)
- Soft headstrap included in delivery

9935-1018737



Laser Safety Glasses (infrared laser light)

Laser Safety glasses for infrared laser light.

Optical Density OD 4 at 780 nm-1100 nm

Protection class DIR L4 at 780 nm-1100 nm

Optical Density OD 5 at 795 nm-1085 nm

Protection class DIR L5 at 795 nm-185 nm

9935-1004887




Therapist Safety Glasses for Red Light

Laser safety glasses for red laser light (660nm)

9935-1008815

396 acupuncture points





Easy and accurate identification of acupuncture points

Conduction paths are coloured according to the 5 phases

3B Male & Female Acupuncture Models

These acupuncture models set the standard in didactics and aesthetics and have been developed in co-operation with experienced European acupuncture specialists. They are about 70 cm tall and contain 396 acupuncture points. All conduction paths are coloured according to the 5 phases. We have intentionally dispensed with Chinese labelling on the model. This makes the accurate identification of acupuncture points significantly easier. True to the motto: less is sometimes more.


-  D/E/F//S Male 9935-1000378
-  D/E/F//S Female 9935-1000379



3B Acupuncture Ears with Needles

These slightly larger than life-size replicas of the human auricula (ear concha) facilitate learning and practising ear acupuncture (also auricular medicine or auricular therapy). The models are made of high-quality 3B SKINlike™ silicone – your guarantee for long life and a realistic ear surface structure. 18 anatomical structures are cited in an enclosed product manual to aid in finding the acupuncture points.

Ear, left	9935-1000374
Ear, right	9935-1000375
Ear set with needles	9935-1000373
Set of 10, left and right	9935-1000376



Ear Acupuncture – 50 x 67 cm

Paper	9935-4006731
Laminated	9935-1001628



Map of Ear Acupuncture

Points from the Chinese and French School of Ear acupuncture, Nogier Frequencyzones, Tables of Frequencies (Nogier, Bahr, Reininger). English. Size DIN A1.

Paper	9935-1013827
-------	--------------


High aluminium rod (top and bottom), printed on heavy, high quality paper.

Canvas	9935-1015542
--------	--------------



Foot reflexology – 50 x 67 cm

Paper	9935-4006729
Laminated	9935-1001624



Body Acupuncture Chart – 98 x 68 cm

Paper	9935-4006730
Laminated	9935-1001626



3BTAPE



3BTAPE is made from a cotton rayon fibre that is **latex free and is hypoallergenic**. **3BTAPES** have a skin-friendly, thermoplastic wavelike acrylic adhesive on one side. The wave like adhesive is used to push sweat and moisture away from the tape. This means that **3BTAPES** are water resistant, air-permeable and do not come off in the shower, pool, or sauna. They can be stretched lengthways by approximately 30 - 40%.

3BTAPE rolls are 5m x 5cm in length and are available in the following colours:

Colour	Item No.
Black	9935-1008621
Beige	9935-1008620
Blue	9935-1008622
Pink	9935-1013842
Yellow	9935-1012803
Green	9935-1012804



New **3BTAPE** Bulk rolls are perfect for clinics and athletic trainers who are constantly using kinesiology tape on their patients who want to buy in bulk and save money. 31 m x 5cm

Colour	Item No.
Black	9935-1013840
Beige	9935-1013836
Blue	9935-1013841
Pink	9935-1013842

Save time with ProCuts!

3BTAPE ProCut Strips

When using kinesiology taping in sports or therapeutic treatment there are four commonly used shapes:



I Strips – 40 p/pack

Beige	9935-1014600
Black	9935-1014602
Blue	9935-1014601
Pink	9935-1014603



X Strips – 40 p/pack

Beige	9935-1014604
Black	9935-1014607
Blue	9935-1014605
Pink	9935-1014608



Y Strips – 40 p/pack

Beige	9935-1014609
Black	9935-1014612
Blue	9935-1014611
Pink	9935-1014613



Eyelash Strips – 40 p/pack

Beige	9935-1014593
Black	9935-1014595
Blue	9935-1014596
Pink	9935-1014598



3BTAPE Coated Kinesiology Scisso

The latest innovation in engineering and craftsmanship, the **3BTAPE Coated Kinesiology Taping Scissors**. **3BTAPE Coated Pro Kinesiology Taping Scissors** are made from high carbon stainless steel which keeps the blade razor-sharp for years. These Kinesiology Taping Scissors are specially coated with black carbon and a fluorine resin which prevents the scissors from getting gummed up by the tape adhesive. Easily cuts multiple strips of kinesiology tape at the same time. Made from premium high carbon stainless steel, the blade stays razor sharp for years. Special coating virtually eliminates adhesive residue. Perfectly balanced for precision cutting.

9935-1014148



Spiral Power Tapes

The spiral tapes offered here are fixed to painful points, scars, trigger points, and acupuncture points. Spiral tapes enable energy to flow away from blockages more effectively, and ensure that the energy flows smoothly. After even short exposure times, many notice an improvement in their mobility and a reduction in pain.

The material is made of cotton fabric, has a thermoplastic acrylic adhesive, and is latex-free. Our spiral tapes are well-tolerated by the skin, are water-resistant, breathable and are easy to remove. They can be worn for several days at a time.

Spiral Power Tapes come in packs of: 180 (22 x 27 mm), 120 (28 x 36 mm), and 40 (28 x 36 mm).

Pack of 180	9935-1013277
Pack of 120	9935-1013276
Pack of 40	9935-1013275




IonFlow 50 Air Purifier

IonFlow air purifier is ideal for medical practices, offices, bedrooms, children's rooms and living rooms. What is so special about IonFlow?

- Effective removal of even the smallest harmful particles from the air in the room
- Greater feeling of wellbeing thanks to healthy air
- Ozone-free (not measurable)
- Very quiet (under 5 dB)
- Energy-saving (under 7 W)
- No need for a filter change
- In just 5 hours, more than 99% of the smallest and most harmful substances have been removed

A. Surface	9935-1015674
B. Style	9935-1015675
C. Signature	9935-1015676



Set up your treatment table anywhere!

Benefits of the portable massage table:

- Maximum working weight: 550 lbs. (250 kg)
- Stable height adjustment from 23 - 33" in. (58 to 85 cm)
- Rubber Feet on legs to prevent scratching and sliding on floor
- Crafted in sturdy natural wood (high quality Birch)



Basic Portable Massage Table

The 3B Basic Portable Massage Table is an economical portable massage table. The table comes standard with 2" dense foam upholstered in a very comfortable vinyl that is oil and waterproof. Includes adjustable face cradle base and pad, side arm extenders (adds 4" to each side of table), hanging arm hammock and standard one pocket carrying case.

Blue	9935-1013724
Green	9935-1013725
Burgundy	9935-1013726



Deluxe Portable Massage Table

The 3B Deluxe Portable Massage Table is a high quality massage table. This table comes with 3" high density small foam and is covered in a comfortable PU vinyl. Includes adjustable face cradle base and pad, hanging arm hammock and standard one pocket carrying case.

Blue	9935-1013727
Green	9935-1013728
Burgundy	9935-1013729

Optional equipment for the portable massage tables, choose your colour:




Adjustable Headrest with Metal Brackets

Blue	9935-1013732
Green	9935-1013734
Burgundy	9935-1013733



Arm Support

Blue	9935-1013738
Green	9935-1013740
Burgundy	9935-1013739



Armrests with Wood Bracket

Blue	9935-1013735
Green	9935-1013737
Burgundy	9935-1013736

MARCA LAB

Alverca Park - Núcleo Empresarial
Armazém Fracção A · Estrada Nacional 10, Km 126,5
2615-701 Sobralinho - Alverca do Ribatejo · Portugal
Tlf.: +351 219 598 163
www.marcalab.pt · geral@marcalab.pt

SEIRIN[®]

Use SEIRIN[®] needles
for gentle, painless,
and stress free

ACUPUNCTURE!

